



6143 North 60th Street
Milwaukee, WI 53218
www.chimerx.com

Certificate of Analysis

Telephone: (414) 535-9506
Toll Free: (800) 626-7833
Fax: (414) 535-9508
Web Support: chimerx.com/support
Email: info@chimerx.com

Bgl II

Catalog No: 2090

Lot No: 2404024

Package Size: See Product Label

Concentration: 10 units/ μ l

Source: *Bacillus globigii*

Storage Conditions: Store at -20°C

Recognition/Cut Site

5'-GCCNNNN↓NGGC-3'

3'-CGGN↑NNNNCCG-5'

Unit Definition

One unit is the amount of enzyme required to completely digest 1 μ g of Lambda DNA in 1 hour at 37°C.

Heat Inactivation

65°C for 10 minutes

Reaction Temperature

N/A

Reaction Buffer

50 mM Tris-HCl (pH 7.5 at 37°C)

100 mM NaCl

10 mM MgCl₂

1 mM dithiothreitol (DTT)

100 μ g/ml bovine serum albumin*

*supplied separately

Assay Conditions

50 mM Tris-HCl (pH 7.5 at 37°C)

100 mM NaCl

10 mM MgCl₂

1 mM dithiothreitol (DTT)

100 μ g/ml bovine serum albumin

1 μ g Lambda DNA

Incubation time is 1 hour at 37°C

Reaction volume is 50 μ l

Storage Buffer

50 mM Tris-HCl (pH 7.5)

50 mM KCl

0.1 mM EDTA

1 mM dithiothreitol (DTT)

200 μ g/ml bovine serum albumin

0.1% Triton X-100™

50% (v/v) glycerol

Quality Control

Non-Specific Endonuclease: 20, 35, and 50 units of Bgl II were incubated with Lambda DNA for 5 hours at 37°C and compared to a standard 4 Unit, 1 hour assay.

3'-Exonuclease: 5, 10, 15, 20, and 30 units of Bgl II were incubated with [³H] labeled 3'-ends of Lambda/Taq I DNA fragments for 1 hour at 37°C.

5'-Exonuclease/5'-Phosphatase: 5, 10, 15, 20, and 30 units of Bgl II were incubated with [γ - ³²P] labeled 5' - ends of Lambda/Hae III DNA fragments for 1 hour at 37°C.

Nicking: 20, 35, and 50 units of Bgl II were incubated with Φ X-174 DNA for 5 hours at 37°C.

Ligation/Recut: 4 units of Bgl II per μ g of Lambda DNA were incubated for 1 hour at 37°C. Resulting fragments were ligated with T4 DNA Ligase and recut with Bgl II.

For more information or to re-order, please visit www.chimerx.com